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Drala® Polycarboxylatzement yellowish white

Dental water-based cement Zinc polycarboxylate cement Cementation and cavity lining cement

Preparation

Dosage of powder:

Press measuring stick firmly down into powder. When withdrawn from container, scrape excess powder from stick with spatula. Release the powder onto the glass slab by tapping the measuring stick by means of a spatula.

Dosage of liquid:

Hold the bottle in vertical position and squeeze. Release pressure when the drop separates from the nozzle. Thus small uniform drops are obtained.

Mixing ratio

For cementations: 1 dose of powder + 3 drops of liquid

The recommended by-weight ratio is 1,5 parts of powder to 1 part of liquid.

For linings: 1 dose of powder + 2 drops of liquid

The recommended by-weight ratio is 2,5 - 3,5 parts of powder to 1 part of liquid.

Mixing

In order to achieve correct mixing of powder and liquid, a glass slab with a temperature of 18 °C - 21 °C (64 °F - 70 °F) and a stainless steel spatula should be used.

Be careful to use only clean glass slab and spatula!

Use spatulation technique that involves placing the powder into the liquid at one time. To facilitate adaptation of a thick and firm mixing, the instrument may be dipped in Drala® Polycarboxylate Cement powder or in alcohol.

General considerations

- The prepared polycarboxylate cement appears thicker than other cements, its flow properties remain excellent however.
- If dosage is made visually, care must be taken to add sufficient powder. Too small amounts of powder retard setting time and may impair the quality of the cement.
- Fresh amounts of liquid never should be added during mixing!
- Drala® Polycarboxylate cement liquid should be dispensed immediately before mixing. Exposure to air leads to loss of water by evaporation.
- Drala® Polycarboxylate cement adheres to tooth structures and stainless steel only if applied while the mix is still glossy. Don't use the mixture when it has lost its gloss.
- The best adhesion is achieved by applying to dry tooth structure. For maximum adhesion tooth surfaces and items to be cemented must be clean and free from contamination. After trial in mouth remove saliva carefully again.
- If possible, remove excess cement before its setting. In order to facilitate the removal of hardened Drala® Polycarboxylate cement excess, the areas adjacent to the cementation could be insulated with some petrolatum.
- It is advisable to clean the instruments with water before the material has set. After complete setting Drala® Polycarboxylate cement can be removed by means of 10 % solution of sodium hydroxide.

Storage

Close recipients carefully after each use. Keep tightly closed at temperatures of 10 °C - 24 °C/50 °F - 75 °F.

Keep away from direct sunlight.

Indications

- Fixing cement for crowns, bridges, inlays.
- Cavity lining cement to protect against thermal and chemical irritants, if a sufficient pulp protection is existing.

Technical Data

DIN EN ISO 9917-1

Contra indications

- The cement is not allowed to be applied directly in close proximity to the unprotected pulp.

Attention

In case of deep cavities a calcium hydroxide lining should be applied to the cavity before applying the zinc polycarboxylate cement.

Mixing time

- 30 seconds

Working time

- **For cementations:** 2 ½ minutes*
- **For linings:** 1 ½ minutes*

Setting time

- **For cementations:** 7 - 8 minutes**
- **For linings:** 5 - 6 minutes**

* from beginning of mixing at 21 °C/70 °F; 50 ± 5% rel. humidity

** from end of mixing at 37 °C/99 °F; 50 ± 5% rel. humidity

Please note

Overdosing of powder or increased temperatures accelerate the setting time.

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GmbH & Co. KG